


Setup for Genotyping Assays:

- Open the **LightCycler® 480 software** and login with your username and password.
 - Choose **New Experiment**.
- Define your PCR program in the **Run Protocol** tab.
 - Select **Dual Color Hydrolysis Probe / UPL Probe** as Detection Format.
 - Select a reaction volume of **20 µl** and setup the **Program**:

Program Name	Cycles	Analysis Mode	Target (°C)	Acquisition Mode	Hold (hh:mm:ss)	Ramp Rate (°C/s)
pre-incubation	1	None	95°C	none	00:03:00	4.4
amplification	40	Quantification	95°C	none	00:00:15	4.4
			60°C	single	00:01:00	2.2
cooling	1	None	40°C	none	00:00:30	2.2

- Add program steps in the **Programs** window with "+" and edit **Cycles** and **Analysis Mode**. To edit **Target (°C)**, **Acquisition Mode** and **Hold**, click on the corresponding step in the **Program Name** window and change parameters or add steps ("+" in the **Temperature Targets** window.
- Click on **Subset Editor** button on the left side of the window.
 - Press "+" to create a new subset and rename your subset.
 - Select wells in the grid and press **Apply**.
- Click on **Sample Editor** button on the left side of the window.
 - **Select Workflow > Endpt Geno**.
 - **Select Filter Combinations > 465-510 nm (FAM) and 533-580 nm (HEX)**.
 - Choose your **Subset** of Samples.



- Define your **No Template Control (NTC)**:

- **Select Samples** field: select well by mouse-click or two wells by ctrl+mouse click.
- **Edit Endpt Geno Properties** field: Type **NTC** in the **Sample Name** field and press **Enter**. Choose **Negative Control** as **Sample Type**.

Pos	Color	Repl Of	Sample Name	EndPt Sample Type	EndPt Genotype
A1	■		NTC	Negative Con	
A2	■		NTC	Negative C ▾	

- Define your **Positive Controls** as Standard (*alternatively: Define your Positive Controls as Unknown*):

- **Select Samples** field: select a well by mouse click.
- **Edit Endpt Geno Properties** field: Type the name of your **Positive Control** for **HEX** in the **Sample Name** field. Choose **Standard**.
- In the chart (EndPt Genotype) type in the genotype corresponding to your **Positive Control HEX**. This will be in most cases **wild type (WT)** for the fluorophore HEX.

Pos	Color	Repl Of	Sample Name	EndPt Sample Type	EndPt Genotype
B3	■		Pos Ctrl HEX	Standard	WT
B4	■		Pos Ctrl FAM	Standard	MUT ▾

- Repeat the steps above with your **Positive Control FAM**. The genotype will be in most cases **mutant (MUT)** for the fluorophore FAM.


- Define your **Samples**:


- **Select Samples** field: select a well by mouse click.
- **Edit Endpt Geno Properties** field: Type the name of your first sample in the corresponding field and press **Enter**. Check **Unknown**.
- Define the rest of your samples.

- **Save** or **export** your experiment by pressing the corresponding button:
- Load your samples and start the experiment.





Analysis of Genotyping Assays:

- Open the **LightCycler® 480 software** and login with your username and password. The **Overview** window appears.
- Click on  or choose **Navigator** in the flip-window on the top left.
 - Choose an experiment from the data bank, or
 - import an experiment located outside the data bank by pressing **Import**.
- After the file is loaded the **Summary** window of your experiment is displayed.
- Press the **Analysis** button to reach the analysis window.
 - Within the **Create New Analysis** field choose **Endpoint Genotyping**.
 - A pop-up window will be launched. If applicable select a **Subset** of samples or analyze **All Samples** in case your plate contains only one type of assay. Give a **Name** to your analysis. Press the **OK** button.



Create new analysis	
Analysis Type *	Endpoint Genotyping
Subset *	MTHFR-A1298C
Program *	amplification
Name *	Endpoint Genotyping for MTHFR-A1298C

- Assign **Allele X** to **FAM** (465-510) and **Allele Y** to **HEX** (533-580) and press **OK**.
 - A Scatter Plot with the signals for HEX (y-axis) and FAM (x-axis) is displayed.
*Optional: Press **Color Comp** and choose **In Use** or – if available – **In Database** for color compensation of FAM (510) and VIC (580). Press the **OK** button.*
 - Press the **Calculate** button and review your results in the chart bottom-left. The column **Results > Call** displays the genotype of your samples.
In the **Scatter Plot** points for homozygous samples for the allele Y (HEX) group along the vertical axis, homozygous samples for allele X (FAM) along the horizontal axis. Heterozygotes will generate a cluster in the middle.
- **Save** or **Export** your data by pressing the corresponding button.  or 
 - After saving your data you can customize and generate a report via the **Report** button.